

## Agilent xCELLigence RTCA DP System

For label-free, real-time cellular analysis

## Introduction

The Agilent xCELLigence real-time cell analysis (RTCA) dual purpose (DP) system provides a unique and powerful means to monitor cells in real time without the potential artifacts generated by using labels. This noninvasive measurement allows detection of changes in adherence, morphology, and viability without needing to overexpress reporter and target proteins. This provides highly physiologically relevant data throughout the experiment.

The Agilent xCELLigence E-Plate features an innovative biosensor configuration that covers 80% of each well bottom surface area. The real-time measurement of impedance across the biosensors provides sensitive immediate detection of the cellular condition and response. This enables a wide array of potential applications including (but not limited to):

- · Cell proliferation
- Cell quality
- Compound-mediated cytotoxicity
- · Cell-mediated cytotoxicity
- Cell adhesion and spreading
- Functional monitoring of receptor tyrosine kinase and GPCR signaling
- · Cell-mediated cytolysis
- Barrier function
- Viral quantification

The Agilent xCELLigence CIM-Plate is a modified Boyden chamber design comprised of a disposable top and bottom chamber featuring the same innovative microelectrode configuration on the bottom of a microporous polyethylene terephthalate (PET) membrane. The median pore size of this membrane is 8  $\mu m$ . The design and real-time measurement allow for precise, quantitative, and dynamic information for applications in cell invasion and migration.

E-Plate 16		
Dimensions	4.0 cm × 8.7 cm × 1.96 cm (W × D × H, with cover)	
Well spacing	9 mm center-to-center as per ANSI/SBS 4-2004 standard	
Well volume	270 ±10 μL	
Well bottom diameter	5.0 ±0.075 mm	
Electrical interface	Compatible with RTCA DP analyzer	
Sensor impedance	$17\pm 5\Omega$ at 10 kHz, when measured with a 1x PBS solution	
Materials	Polystyrene well plate, glass sensor substrate, UV irradiated	
Environment	Temperature: +15 to +40 °C, relative humidity: 98% maximum without condensation	

E-Plate VIEW 16		
All E-Plate 16 specifications apply		
Viewing window	Four center electrodes removed to aid in microscopy (~400 µm width)	

CIM-Plate 16		
Dimensions	4.0 cm × 8.7 cm × 2.6 cm (W × D × H, assembled, with cover)	
Well spacing	9 mm center-to-center as per ANSI/SBS 4-2004 standard	
Upper well volume	180 ±5 μL	
Lower well volume	162 ±3 μL	
Membrane	PET membrane with 8 µm pore size	
Well bottom diameter	5.0 ±0.075 mm	
Electrical interface	Compatible with RTCA DP analyzer	
Sensor impedance	24 ±8 Ω at 10 kHz, when measured with a 1x PBS solution	
Materials	PET well plate, PET membrane sensor substrate UV irradiated	
Environment	Temperature: +15 to +40 °C, relative humidity: 98% maximum without condensation	



RTCA DP Analyzer		
Electrical input	+5 VDC, 1 W max.	
Electronic switch resistance	2 to 5 Ω	
Electronic interface	Handling three E-Plate 16 devices or three CIM-Plate devices	
Communication	USB 2.0	
Environment	Temperature: +20 to +40 °C, relative humidity: 98% max noncondensing	
Output test signal	22 mV rms ±(2% +5 mVrms) at 10, 25, and 50 kHz	
Impedance measure accuracy	±(1% + 1.5 Ω)	
Impedance measurement repeatability	0.8%	
Impedance dynamic range	10 to 5,000 Ω	
Status indicators	Analyzer status	



RTCA DP Control Unit	
Computer with pre-installed RTCA software	
User-friendly graphical user interface (GUI)	
≥500 GB hard disk drive	
≥4 GB RAM	
≥256 MB graphics device	

## www.agilent.com/chem

For Research Use Only. Not for use in diagnostic procedures.

This information is subject to change without notice.

© Agilent Technologies, Inc. 2019 Printed in the USA, November 1, 2019 5994-1068EN DE.5735648148

